

Anti-Ranbp1 Polyclonal Antibody For Wb, Elisa - P43487

Item Number: CM0020546



Introduction

Rabbit polyclonal antibody targeting human RANBP1 (RanBP1/HTF9A) for Western Blot and ELISA. Validated in human, mouse, and rat. Ideal for nucleocytoplasmic transport and mitosis research.

[Learn More](#)

Parameter	Value
Product Name	RANBP1 Rabbit pAb
Remarks/Alias	HTF9A; RANBP1
Species	Human
GeneID (Human)	5902
GeneID	5902
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human RANBP1 (NP_002873.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P43487
Protein Weight	23kDa
Shipping	Ice bag