

Anti-Ranbp10 Rabbit Polyclonal Antibody For Wb And Elisa - Q6Vn20

Item Number: CM0019114

Introduction

Rabbit polyclonal antibody targeting human RANBP10, validated for Western blot and ELISA. Cross-reacts with human and mouse. Suitable for studying CTLH E3 ubiquitin ligase function, Ran GTPase signaling, and platelet biology.

[Learn More](#)

Parameter	Value
Product Name	RANBP10 Rabbit pAb
Remarks/Alias	RANBP10
Species	Human
GeneID (Human)	57610
GeneID	57610
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 362-497 of human RANBP10 (NP_065901.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q6VN20
Protein Weight	67kDa
Shipping	Ice bag