


Anti-Ranbp2 Polyclonal Antibody For Wb, Ihc-P, Elisa - P49792

Item Number: CM0023894



Introduction

Rabbit polyclonal antibody targeting human RANBP2 (NUP358), validated for WB, IHC-P, ELISA; cross-reacts with mouse and rat. Ideal for nucleocytoplasmic transport and SUMOylation research.

[Learn More](#)

Parameter	Value
Product Name	RANBP2 Rabbit pAb
Remarks/Alias	ANE1; TRP1; TRP2; ADANE; IIAE3; NUP358; RANBP2
Species	Human
Gene ID (Human)	5903
Gene ID	5903
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 2900-3050 of human RANBP2 (NP_006258.3))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P49792
Protein Weight	358kDa
Shipping	Ice bag