

Anti-Ranbp9 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q96S59

Item Number: CM0003791



Introduction

Rabbit monoclonal antibody targeting human RanBP9 (BPM-L/RanBPM). Validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Suitable for research on ubiquitination, signaling, and cell adhesion.

[Learn More](#)

Parameter	Value
Product Name	RanBP9 Rabbit mAb
Remarks/Alias	BPM-L; BPM90; RANBPM; RanBP7; RanBP9
Species	Human
Gene ID (Human)	10048
Gene ID	10048
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human RanBP9 (Q96S59)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96S59
Protein Weight	78kDa
Shipping	Ice bag