

Anti-Ranbp9 Rabbit Polyclonal Antibody For Wb, Elisa - Q96S59

Item Number: CM0018394



Introduction

High-quality rabbit polyclonal antibody targeting human RanBP9 (RanBPM). Validated for WB and ELISA, this antibody is ideal for studying CTLH E3 ligase complex function, cell signaling, and cytoskeletal rearrangement. Reacts with human samples.

[Learn More](#)

Parameter	Value
Product Name	RanBP9 Rabbit pAb
Remarks/Alias	BPM-L; BPM90; RANBPM; RanBP7; RanBP9
Species	Human
Gene ID (Human)	10048
Gene ID	10048
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 480-605 of human RanBP9 (NP_005484.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q96S59
Protein Weight	78kDa
Shipping	Ice bag