

Anti-Rangap1 Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P46060

Item Number: CM0011379



Introduction

Recombinant rabbit monoclonal antibody against human RanGAP1 (P46060), validated for WB, IHC-P, IF/ICC, and ELISA. Recognizes mouse and rat. For Ran GTPase activation and nuclear transport studies.

[Learn More](#)

Parameter	Value
Product Name	RanGAP1 Rabbit mAb
Remarks/Alias	SD; Fug1; RANGAP; RanGAP1
Species	Human
GeneID (human)	5905
GeneID	5905
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 1-100 of human RanGAP1 (P46060).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P46060
Protein Weight	64kDa
Shipping	Ice bag