


Anti-Rankl Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O14788

Item Number: CM0017713



Introduction

High-quality rabbit polyclonal antibody specific to human RANKL (TNFSF11), validated for WB, IF/ICC, and ELISA applications. Cross-reacts with mouse and rat. Ideal for osteoclast differentiation and immune response studies. - O14788

[Learn More](#)

Parameter	Specification
Product Name	RANKL Rabbit pAb
Remarks/Alias	ODF; OPGL; sOdf; CD254; OPTB2; RANKL; TNLG6B; TRANCE; hRANKL2
Species	Human
GeneID (Human)	8600
GeneID	8600
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 78-317 of human RANKL (NP_003692.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14788
Protein Weight	35kDa
Shipping	Ice bag