

Anti-Rap1A Polyclonal Antibody For Wb, Ihc-P, Elisa - P62834

Item Number: CM0016760



Introduction

KO-validated anti-Rap1A rabbit polyclonal antibody detects human, mouse, rat Rap1A in WB, IHC-P, ELISA. Ideal for studies on Ras counteraction, neurite outgrowth, and vascular biology.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Rap1A Rabbit pAb
Remarks/Alias	RAP1; C21KG; G-22K; KREV1; KREV-1; SMGP21; 1A
Species	Human
GeneID (human)	5906
GeneID	5906
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 39-138 of human Rap1A (NP_002875.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62834
Protein Weight	21kDa
Shipping	Ice bag