

Anti-Rap1A + Rap1B Monoclonal Antibody For Wb, If/Icc, Elisa - P61224 / P62834

Item Number: CM0012738



Introduction

Recombinant rabbit monoclonal antibody targeting human RAP1A and RAP1B, validated for Western blot, IF/ICC, and ELISA. Ideal for studying endothelial cell polarity and vascular lumen formation.

[Learn More](#)

Parameter	Value
Product Name	RAP1A + RAP1B Rabbit mAb
Remarks/Alias	K-REV, RAL1B, RAP1A + RAP1B
Species	Human
Gene ID (Human)	5906 (RAP1A) / 5908 (RAP1B)
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 70-150 of human RAP1A + RAP1B (P61224)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61224 / P62834
Protein Weight	21kDa
Shipping	Ice bag