

Anti-Rap1B Polyclonal Antibody For Wb, If/Icc, Elisa - P61224

Item Number: CM0016868



Introduction

Rabbit polyclonal antibody targeting human Rap1B (P61224) for WB, IF/ICC, ELISA. Validated in human and mouse samples. Ideal for endothelial cell biology and GTPase signaling research.

[Learn More](#)

Parameter	Specification
Product Name	Rap1B Rabbit pAb
Remarks/Alias	K-REV; RAL1B; Rap1B
Species	Human
Gene ID (Human)	5908
Gene ID	5908
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-184 of human Rap1B (NP_056461.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P61224
Protein Weight	21kDa
Shipping	Ice bag