

Anti-Rap1B Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P61224

Item Number: CM0019954



Introduction

Rabbit polyclonal antibody targeting RAP1B (Ras-related protein Rap-1b). Validated for WB, IF/ICC, ELISA; reacts with human, mouse, rat. Ideal for studying endothelial cell polarity, vascular biology, and GTPase signaling.

[Learn More](#)

Parameter	Value
Product Name	RAP1B Rabbit pAb
Remarks/Alias	RAP1B; K-REV; RAL1B; ras-related protein Rap-1b
Species	Human
Gene ID (Human)	5908
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 81-181 of human RAP1B (NP_056461.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61224
Protein Weight	15kDa/16kDa/18kDa/20kDa
Shipping	Ice bag