


Anti-Rapgef1 (C3G/Grf2) Polyclonal Antibody For Wb, If/Icc, Elisa - Q13905

Item Number: CM0022436



Introduction

Rabbit polyclonal antibody against human RAPGEF1/C3G (Q13905), validated for WB, IF/ICC, ELISA. Ideal for studying cell adhesion, RAS activation, and neurite outgrowth in human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	RAPGEF1 Rabbit pAb
Remarks/Alias	C3G; GRF2; RAPGEF1
Species	Human
GeneID (Human)	2889
GeneID	2889
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 844-1095 of human RAPGEF1 (NP_941372.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13905
Protein Weight	121kDa
Shipping	Ice bag