

Anti-Raph1 Rabbit Pab For Wb, If/Icc, Elisa - Q70E73

Item Number: CM0015576



Introduction

Anti-RAPH1 polyclonal antibody validated for WB, IF/ICC, and ELISA in human, mouse, and rat. Explore lamellipodial dynamics and cell adhesion mechanisms with this IVD raw material from CamelBio.

[Learn More](#)

Parameter	Value
Product Name	RAPH1 Rabbit pAb
Remarks/Alias	LPD; RMO1; PREL2; PREL-2; ALS2CR9; ALS2CR18; RaIGDS/AF-6; RAPH1
Species	Human
GeneID (Human)	65059
GeneID	65059
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human RAPH1 (NP_998754.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q70E73
Protein Weight	135kDa
Shipping	Ice bag