

Anti-Rapsn Polyclonal Antibody For Wb, Ihc-P, Elisa - Q13702

Item Number: CM0028661

Introduction

High-quality rabbit polyclonal antibody against human RAPSN, validated for Western blot, IHC-P, and ELISA. Detects endogenous levels of RAPSN in human, mouse, and rat samples. Ideal for neuromuscular junction research.

[Learn More](#)

Parameter	Value
Product Name	RAPSN Rabbit pAb
Remarks/Alias	FADS; CMS11; CMS4C; FADS2; RAPSIN; RNF205; RAPSN
Species	Human
Gene ID (Human)	5913
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-353 of human RAPSN (NP_116034.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13702
Protein Weight	46kDa
Shipping	Ice bag