

Anti-Raptor Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q8N122

Item Number: CM0012729



Introduction

CamelBio's Raptor Rabbit mAb targets human/mouse Raptor (Q8N122) for WB, IHC-P, ELISA. This mTORC1 scaffold protein is key to nutrient sensing, autophagy, and cell growth signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	Raptor Rabbit mAb
Remarks/Alias	KOG1; Mip1; Raptor
Species	Human
Gene ID (Human)	57521
Gene ID	57521
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1236-1335 of human Raptor (Q8N122).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8N122
Protein Weight	149kDa
Shipping	Ice bag