

Anti-Raptor Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q8N122

Item Number: CM0023567



Introduction

Rabbit polyclonal antibody targeting Raptor (RPTOR), a key mTORC1 component. Validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Suitable for cell growth, nutrient signaling, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Raptor Rabbit pAb
Remarks/Alias	KOG1; Mip1; Raptor
Species	Human
Gene ID (Human)	57521
Gene ID	57521
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 888-1013 of human RPTOR (NP_065812.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N122
Protein Weight	149kDa
Shipping	Ice bag