


Anti-Rarg Rabbit Monoclonal Antibody For Wb, Elisa - P13631

Item Number: CM0004906



Introduction

Rabbit monoclonal antibody against human RARG, validated for Western blot and ELISA. Cross-reacts with human, mouse, and rat RARG. Suitable for studying retinoic acid signaling and gene regulation.

[Learn More](#)

Parameter	Value
Product Name	RARG Rabbit mAb
Remarks/Alias	RARC; NR1B3; RARGgamma; RARG
Species	Human
Gene ID	5916
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human RARG (NP_000957.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13631
Protein Weight	50kDa
Shipping	Ice bag