


Anti-Rarg Rabbit Polyclonal Antibody For Wb, Elisa - P13631

Item Number: CM0030722



Introduction

High-quality anti-RARG rabbit polyclonal antibody validated for WB and ELISA. Recognizes human, mouse, and rat RARgamma. Ideal for studying retinoic acid signaling, gene regulation, and developmental biology.

[Learn More](#)

Parameter	Value
Product Name	RARG Rabbit pAb
Remarks/Alias	RARC; NR1B3; RARgamma; RARG
Species	Human
GeneID (Human)	5916
GeneID	5916
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 168-199 of human RARG (NP_000957.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13631
Protein Weight	50kDa
Shipping	Ice bag