

Anti-Rarres2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q99969

Item Number: CM0030954



Introduction

Rabbit polyclonal antibody targeting human RARRES2 (chemerin/TIG2) for Western blot, IF/ICC, and ELISA. Cross-reacts with mouse. Suitable for adipokine research in metabolism, inflammation, and obesity studies.

[Learn More](#)

Parameter	Value
Product Name	RARRES2 Rabbit pAb
Remarks/Alias	TIG2; HP10433; RARRES2
Species	Human
GeneID (Human)	5919
GeneID	5919
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 17-163 of human RARRES2 (NP_002880.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q99969
Protein Weight	19kDa
Shipping	Ice bag