


Anti-Rars Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P54136

Item Number: CM0009724



Introduction

Rabbit monoclonal antibody against human RARS (arginyl-tRNA synthetase), validated for WB, IF-P, IHC-P, ELISA. Cross-reacts with mouse. Ideal for protein synthesis and inflammatory cytokine studies.

[Learn More](#)

Parameter	Value
Product Name	RARS Rabbit mAb
Remarks/Alias	HLD9; RARS; ArgRS; DALRD1
Species	Human
Gene ID (Human)	5917
Gene ID	5917
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human RARS (NP_002878.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P54136
Protein Weight	75kDa
Shipping	Ice bag