


Anti-Rars Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P54136

Item Number: CM0031035



Introduction

Rabbit polyclonal antibody against human RARS (ArgRS). Validated for WB, IHC-P, IF/ICC, ELISA, with cross-reactivity to human, mouse, and rat. Ideal for protein synthesis and aminoacyl-tRNA synthetase research.

[Learn More](#)

Parameter	Value
Product Name	RARS Rabbit pAb
Remarks/Alias	HLD9; RARS; ArgRS; DALRD1
Species	Human
GeneID (human)	5917
GeneID	5917
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human RARS (NP_002878.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P54136
Protein Weight	75kDa
Shipping	Ice bag