

Anti-Rarβ Monoclonal Antibody For Wb, If/Icc, Elisa - P10826

Item Number: CM0005138



Introduction

Rabbit monoclonal antibody against human RARβ (retinoic acid receptor beta). Validated for WB, IF/ICC, and ELISA. Cross-reacts with human and mouse. Ideal for studying retinoic acid signaling and nuclear receptor function.

[Learn More](#)

Parameter	Value
Product Name	RARβ Rabbit mAb
Remarks/Alias	HAP; RRB2; NR1B2; MCOP512; RARbeta; RARbeta1; RARβ
Species	Human
GeneID (Human)	5915
GeneID	5915
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 350-455 of human RARβ (NP_000956.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P10826
Protein Weight	50kDa
Shipping	Ice bag