


Anti-Rasa1 Monoclonal Antibody For Wb, Elisa - P20936

Item Number: CM0009035



Introduction

Rabbit monoclonal antibody against human RASA1 (Ras GTPase-activating protein 1, p120GAP) for WB and ELISA. Cross-reacts with mouse. Ideal for Ras signaling research. Target molecular weight ~116 kDa.

[Learn More](#)

Parameter	Value
Product Name	RASA1 Rabbit mAb
Remarks/Alias	GAP, PKWS, RASA, p120, CMAVM, CM-AVM, CMAVM1, RASGAP, p120GAP, p120RASGAP, RASA1
Species	Human
GeneID (Human)	5921
GeneID	5921
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human RASA1 (P20936)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P20936
Protein Weight	116kDa
Shipping	Ice bag