

Anti-Rasal1 Rabbit Polyclonal Antibody For Wb And Elisa - 095294

Item Number: CM0023298

Introduction

Anti-RASAL1 rabbit polyclonal antibody validated for Western blot (WB) and ELISA. Recognizes human, mouse, and rat RASAL1, a RasGAP inhibitory regulator involved in dendrite formation and cAMP signaling.

[Learn More](#)

Parameter	Value
Product Name	RASAL1 Rabbit pAb
Remarks/Alias	RASAL; RASAL1
Species	Human
GenelD (human)	8437
GenelD	8437
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 120-340 of human RASAL1 (NP_004649.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95294
Protein Weight	90kDa
Shipping	Ice bag