

Anti-Rasal2 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Ujf2

Item Number: CM0016708

Introduction

Anti-RASAL2 Rabbit Polyclonal Antibody validated in WB, ELISA, with cross-reactivity to human and rat. Targets RASAL2, an inhibitory regulator of Ras-cyclic AMP pathway. Sourced from rabbit for research.

[Learn More](#)

Parameter	Value
Product Name	RASAL2 Rabbit pAb
Remarks/Alias	NGAP; RASAL2
Species	Human
GenelD (human)	9462
GenelD	9462
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 750-1040 of human RASAL2 (NP_733793.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9UJF2
Protein Weight	129kDa
Shipping	Ice bag