

Anti-Rasd1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y272

Item Number: CM0015898

Introduction

Anti-RASD1 rabbit polyclonal antibody validated for WB and ELISA in human and mouse. Targets the dexamethasone-induced Ras-related protein 1, a small GTPase involved in transcription regulation.

[Learn More](#)

Parameter	Value
Product Name	RASD1 Rabbit pAb
Remarks/Alias	AGS1; DEXRAS1; MGC:26290; RASD1
Species	Human
GeneID (Human)	51655
GeneID	51655
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 132-281 of human RASD1 (NP_057168.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y272
Protein Weight	32kDa
Shipping	Ice bag