


Anti-Rasgrp1 Polyclonal Antibody For Wb, If/Icc, Elisa - 095267

Item Number: CM0018243



Introduction

Rabbit polyclonal antibody against human RASGRP1, suitable for WB, IF/ICC, and ELISA; cross-reacts with mouse and rat. Useful for studying Ras signaling in immune cell development.

[Learn More](#)

Parameter	Value
Product Name	RASGRP1 Rabbit pAb
Remarks/Alias	IMD64; RASGRP; CALDAG-GEFI; CALDAG-GEFII; RASGRP1
Species	Human
GeneID (Human)	10125
GeneID	10125
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 598-797 of human RASGRP1 (NP_005730.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95267
Protein Weight	90kDa
Shipping	Ice bag