

Anti-Rasgrp1 Rabbit Polyclonal Antibody For Wb, Ip, Elisa - Q9Z1S3

Item Number: CM0028649

Introduction

RASGRP1 rabbit polyclonal antibody for WB, IP, ELISA. Detects human and mouse RASGRP1 (~90 kDa). Suitable for Ras signaling, T/B cell development, and oncogenesis research.

[Learn More](#)

Parameter	Value
Product Name	RASGRP1 Rabbit pAb
Remarks/Alias	Rasgrp; calDAG-GEFII
Species	Mouse
GeneID (Human)	10125
GeneID	19419
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 678-795 of mouse RASGRP1 (NP_035376.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Z1S3
Protein Weight	90kDa
Shipping	Ice bag