


Anti-Rasip1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q5U651

Item Number: CM0024210



Introduction

Rabbit polyclonal antibody against RASIP1, a critical regulator of GTPase signaling in endothelial cells. Suitable for Western blot, immunofluorescence, and ELISA. Cross-reacts with human, mouse, and rat. ~103 kDa.

[Learn More](#)

Parameter	Value
Product Name	RASIP1 Rabbit pAb
Remarks/Alias	RAIN; RASIP1
Species	Human
GeneID (Human)	54922
GeneID	54922
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 635-785 of human RASIP1 (Q5U651).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5U651
Protein Weight	103kDa
Shipping	Ice bag