


Anti-Rassf1 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Ns23

Item Number: CM0017055



Introduction

High-quality rabbit polyclonal antibody against human RASSF1 (UniProt Q9NS23), validated for WB and ELISA. Suitable for tumor suppressor research, cell cycle and apoptosis studies. Reacts with human and mouse.

[Learn More](#)

Parameter	Description
Product Name	RASSF1 Rabbit pAb
Remarks/Alias	123F2; RDA32; NORE2A; RASSF1A; REH3P21; RASSF1
Species	Human
GeneID (Human)	11186
GeneID	11186
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 101-340 of human RASSF1 (NP_009113.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NS23
Protein Weight	39kDa
Shipping	Ice bag