

Anti-Rassf4 Polyclonal Antibody For Wb, Elisa - Q9H2L5

Item Number: CM0016607

Introduction

Validate RASSF4, a potential tumor suppressor and KRAS effector, with this rabbit polyclonal antibody. Suitable for WB and ELISA, cross-reacts with human, mouse, and rat. Ideal for apoptosis and cell cycle studies.

[Learn More](#)

Parameter	Value
Product Name	RASSF4 Rabbit pAb
Remarks/Alias	AD037; RASSF4
Species	Human
GeneID (human)	83937
GeneID	83937
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human RASSF4 (NP_114412.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H2L5
Protein Weight	37kDa
Shipping	Ice bag