

Anti-Rassf5 Polyclonal Antibody For Wb, Elisa - Q8Www0

Item Number: CM0031087



Introduction

RASSF5 rabbit polyclonal antibody validated for WB and ELISA, cross-reactive with human, mouse, rat. Targets the tumor suppressor NORE1/RAPL (47 kDa) involved in apoptosis and lymphocyte adhesion.

[Learn More](#)

Parameter	Value
Product Name	RASSF5 Rabbit pAb
Remarks/Alias	RAPL; Maxp1; NORE1; NORE1A; NORE1B; RASSF3; RASSF5
Species	Human
GeneID (Human)	83593
GeneID	83593
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-265 of human RASSF5 (NP_872606.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WWW0
Protein Weight	47kDa
Shipping	Ice bag