

Anti-Rassf5 Rabbit Polyclonal Antibody For Wb And Elisa - Q8Www0

Item Number: CM0015291



Introduction

RASSF5 Rabbit Polyclonal Antibody for WB and ELISA. Validated in human and mouse samples. Targets the potential tumor suppressor and lymphocyte adhesion regulator RASSF5 (NORE1/RAPL).

[Learn More](#)

Parameter	Value
Product Name	RASSF5 Rabbit pAb
Remarks/Alias	RAPL; Maxp1; NORE1; NORE1A; NORE1B; RASSF3; RASSF5
Species	Human
GeneID (Human)	83593
GeneID	83593
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8WWW0
Protein Weight	47kDa
Shipping	Ice bag