

Anti-Rassf6 Rabbit Polyclonal Antibody For Wb, Elisa - Q6Ztq3

Item Number: CM0029062

Introduction

Rabbit polyclonal antibody targeting RASSF6, suitable for Western blot and ELISA, with cross-reactivity to mouse and rat. Used in apoptosis and tumor biology research.

[Learn More](#)

Parameter	Value
Product Name	RASSF6 Rabbit pAb
Remarks/Alias	RASSF6
Species	Human
GeneID (Human)	166824
GeneID	166824
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-303 of human RASSF6 (NP_001257320.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q6ZTQ3
Protein Weight	43kDa
Shipping	Ice bag