

Anti-Rat C-Reactive Protein/Crp Rabbit Polyclonal Antibody For Wb, Elisa - P48199

Item Number: CM0020688

Introduction

Rabbit polyclonal antibody against rat C-reactive protein (CRP), validated for WB and ELISA. Recombinant fusion protein immunogen (aa 20-230). Recognizes rat CRP (Swiss-Prot P48199). Ideal for host defense and inflammation research.

[Learn More](#)

Parameter	Details
Product Name	Rat C-Reactive Protein/CRP Rabbit pAb
Remarks/Alias	Aa1249; Ac1262; Ab1-341; Ab2-196; Ac1-114; Ac2-069; Ba2-693; Rat C-Reactive Protein/CRP
Species	Rat
Gene ID	25419
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-230 of rat C-Reactive Protein/CRP (NP_058792.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P48199
Protein Weight	25kDa
Shipping	Ice bag