

Anti-Rav1 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Zwm9

Item Number: CM0020576



Introduction

Rabbit polyclonal antibody targeting Arabidopsis thaliana RAV1 transcription factor, suitable for Western blot (WB) and ELISA. Detects a 39 kDa nuclear protein involved in ethylene signaling and growth regulation.

[Learn More](#)

Parameter	Value
Product Name	RAV1 Rabbit pAb
Remarks/Alias	AtRAV1; EDF4; ETHYLENE RESPONSE DNA BINDING FACTOR 4; related to ABI3/VP1 1; T6J4.2; T6J4_2; RAV1
Species	Arabidopsis thaliana
Gene ID	837886
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 245-344 of Arabidopsis thaliana RAV1 (NP_172784.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Arabidopsis thaliana
SWISS	Q9ZWM9
Protein Weight	39kDa
Shipping	Ice bag