


## Anti-Rb Monoclonal Antibody For Ihc-P And Elisa - P06400

Item Number: CM0004049



### Introduction

High-quality rabbit monoclonal antibody against human retinoblastoma protein (RB/pRb) validated for IHC-P and ELISA. Ideal for cancer and cell cycle research, diagnostic development, and target detection.

[Learn More](#)

| Parameter        | Value                                                                                                           |
|------------------|-----------------------------------------------------------------------------------------------------------------|
| Product Name     | RB Rabbit mAb                                                                                                   |
| Remarks/Alias    | RB; pRb; OSRC; pp110; p105-Rb; PPP1R130; p110-RB1                                                               |
| Species          | Human                                                                                                           |
| GeneID (Human)   | 5925                                                                                                            |
| GeneID           | 5925                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 701-928 of human RB (NP_000312.2) |
| Source           | Rabbit                                                                                                          |
| Category         | Monoclonal Antibodies                                                                                           |
| Application      | IHC-P, ELISA                                                                                                    |
| Cross Reactivity | Human                                                                                                           |
| SWISS            | P06400                                                                                                          |
| Protein Weight   | 106kDa                                                                                                          |
| Shipping         | Ice bag                                                                                                         |