


Anti-Rb Polyclonal Antibody For Wb, If/Icc, Elisa - P06400

Item Number: CM0021984



Introduction

Human RB (P06400) rabbit polyclonal antibody validated for WB, IF/ICC, and ELISA.

Recombinant immunogen (aa 701-928) detects full-length ~106 kDa retinoblastoma protein.

Ideal for cell cycle, tumor biology, and host-virus interaction studies.

[Learn More](#)

Parameter	Specification
Product Name	RB Rabbit pAb
Remarks/Alias	RB; pRb; OSRC; pp110; p105-Rb; PPP1R130; p110-RB1
Species	Human
Gene ID (Human)	5925
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 701-928 of human RB (NP_000312.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P06400
Protein Weight	106 kDa
Shipping	Ice bag