

Anti-Rbap46/Rbbp7 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q16576

Item Number: CM0006609



Introduction

RbAp46/RBBP7 Rabbit mAb validated for WB, IHC-P, IF/ICC, ELISA; reacts with human, mouse, rat. Targets histone-binding protein RBBP7, a subunit of chromatin remodeling complexes. Ideal for epigenetic and cancer research.

[Learn More](#)

Parameter	Value
Product Name	RbAp46/RBBP7 Rabbit mAb
Remarks/Alias	RbAp46; RbAp46/RBBP7
Species	Human
Gene ID (Human)	5931
Gene ID	5931
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-230 of human RbAp46/RBBP7 (NP_002884.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16576
Protein Weight	48kDa
Shipping	Ice bag