


Anti-Rbap46/Rbbp7 Rabbit Polyclonal Antibody - Q16576

Item Number: CM0017684



Introduction

Rabbit polyclonal antibody against human RbAp46/RBBP7 histone-binding protein. Validated for WB, IHC-P, IF/ICC, ELISA, and ChIP. Cross-reacts with mouse and rat. Suitable for chromatin research.

[Learn More](#)

Parameter	Value
Product Name	RbAp46/RBBP7 Rabbit pAb
Remarks/Alias	RbAp46; RbAp46/RBBP7
Species	Human
GeneID (Human)	5931
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-230 of human RbAp46/RBBP7 (NP_002884.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16576
Protein Weight	48kDa
Shipping	Ice bag