

Anti-Rbbp4 Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q09028

Item Number: CM0005504



Introduction

Rabbit monoclonal antibody against human RBBP4 (NURF55, RBAP48) for Western blot, immunohistochemistry, immunofluorescence, and ELISA. Validated in human, mouse, and rat samples. Research-use IVD component.

[Learn More](#)

Parameter	Value
Product Name	RBBP4 Rabbit mAb
Remarks/Alias	NURF55; RBAP48; lin-53; RBBP4
Species	Human
Gene ID (Human)	5928
Gene ID	5928
Immunogen	Synthetic peptide corresponding to amino acids 1-100 of human RBBP4 (NP_005601.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q09028
Protein Weight	48kDa
Shipping	Ice bag