

Anti-Rbbp5 Rabbit Monoclonal Antibody For Wb, If, Ihc, Elisa - Q15291

Item Number: CM0006522



Introduction

RBBP5 rabbit monoclonal antibody validated for WB, IF, IHC, ELISA in human, mouse, rat. Targets retinoblastoma-binding protein 5, a core component of MLL/SET1 H3K4 methyltransferase complexes.

[Learn More](#)

Parameter	Value
Product Name	RBBP5 Rabbit mAb
Remarks/Alias	RBQ3; SWD1; RBBP5
Species	Human
Gene ID	5929
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 319-538 of human RBBP5 (NP_005048.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15291
Protein Weight	59kDa
Shipping	Ice bag