


Anti-Rbbp6 Polyclonal Antibody For Wb, Elisa - Q7Z6E9

Item Number: CM0023669



Introduction

Rabbit polyclonal antibody targeting human RBBP6 (E3 ubiquitin-protein ligase), validated for WB and ELISA, with cross-reactivity in mouse and rat. Suitable for p53 regulation, DNA replication, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	RBBP6 Rabbit pAb
Remarks/Alias	PACT; MY038; P2P-R; RBQ-1; SNAMA; RBBP6
Species	Human
Gene ID	5930
Immunogen	Recombinant Protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q7Z6E9
Protein Weight	202kDa
Shipping	Ice bag