

Anti-Rbbp8 Rabbit Polyclonal Antibody For Wb - Q99708

Item Number: CM0014671



Introduction

Rabbit polyclonal antibody against human RBBP8 (CtIP), validated for Western blot and ELISA, with cross-reactivity to rat. Ideal for DNA repair, cell cycle, and oncology research.

[Learn More](#)

Parameter	Value
Product Name	RBBP8 Rabbit pAb
Remarks/Alias	RIM; COM1; CTIP; JWDS; SAE2; JAWAD; SCKL2; RBBP8
Species	Human
GeneID (Human)	5932
GeneID	5932
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 648-897 of human RBBP8 (NP_002885.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q99708
Protein Weight	102kDa
Shipping	Ice bag