


Anti-Rbbp9 Polyclonal Antibody For Wb, Ihc-P, Elisa - 075884

Item Number: CM0028014



Introduction

RBBP9 Rabbit pAb is a polyclonal antibody targeting serine hydrolase RBBP9, validated for WB, IHC-P, and ELISA in human, mouse, rat samples. Ideal for cancer and TGF-beta signaling research.

[Learn More](#)

Parameter	Value
Product Name	RBBP9 Rabbit pAb
Remarks/Alias	BOG; RBBP10; RBBP9
Species	Human
GeneID (Human)	10741
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-186 of human RBBP9 (NP_006597.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075884
Protein Weight	21kDa
Shipping	Ice bag