


Anti-Rbfox1 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - Q9Nwb1

Item Number: CM0030467



Introduction

Rabbit polyclonal antibody against human RBFOX1 (also A2BP1, Fox-1 homolog A) for Western blot, IF-P, and ELISA. Confirmed cross-reactivity with human, mouse, and rat; immunogen N-terminal 130 aa.

[Learn More](#)

Parameter	Value
Product Name	RBFOX1 Rabbit pAb
Remarks/Alias	2BP1; FOX1; A2BP1; FOX-1; HRNBP1; RBFOX1
Species	Human
Gene ID (Human)	54715
Gene ID	54715
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human RBFOX1 (NP_665900.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NWB1
Protein Weight	43kDa
Shipping	Ice bag