

Anti-Rbl2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q08999

Item Number: CM0017064



Introduction

RBL2 Rabbit Polyclonal Antibody (p130, Rb2) validated for WB, IF/ICC, ELISA. Recognizes human and mouse RBL2, a key cell cycle regulator and tumor suppressor. Ideal for chromatin and cancer research.

[Learn More](#)

Parameter	Value
Product Name	RBL2 Rabbit pAb
Remarks/Alias	Rb2; P130; BRUWAG; RBL2
Species	Human
GeneID (Human)	5934
GeneID	5934
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 413-613 of human RBL2 (NP_005602.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q08999
Protein Weight	128kDa
Shipping	Ice bag