


Anti-Rbm14 Polyclonal Antibody For Wb, Elisa - Q96Pk6

Item Number: CM0029353



Introduction

Rabbit polyclonal antibody against human RBM14, suitable for WB and ELISA. Cross-reacts with human, mouse, and rat. Useful for studying RNA-binding protein 14 in transcription and centriole regulation.

[Learn More](#)

Parameter	Value
Product Name	RBM14 Rabbit pAb
Remarks/Alias	SIP; COAA; PSP2; SYTIP1; TMEM137; RBM14
Species	Human
Gene ID (Human)	10432
Gene ID	10432
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 110-220 of human RBM14 (NP_006319.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96PK6
Protein Weight	69kDa
Shipping	Ice bag