


Anti-Rbm14 Polyclonal Antibody For Wb, If/Icc, Elisa - Q96Pk6

Item Number: CM0013255



Introduction

RBM14 Rabbit pAb is a polyclonal antibody raised against human RBM14 (aa 110-220). Validated for WB, IF/ICC, ELISA with reactivity in human, mouse, and rat. Ideal for studying transcriptional regulation, centriole biology, and innate immunity.

[Learn More](#)

Parameter	Value
Product Name	RBM14 Rabbit pAb
Remarks/Alias	SIP; COAA; PSP2; SYTIP1; TMEM137; RBM14
Species	Human
Gene ID	10432
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 110-220 of human RBM14 (NP_006319.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96PK6
Protein Weight	69kDa
Shipping	Ice bag