

Anti-Rbm15 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q96T37

Item Number: CM0026934



Introduction

RBM15 Rabbit Polyclonal Antibody detects human and mouse RBM15 (107 kDa), validated for WB, IHC-P, and ELISA. Ideal for studying m6A methylation, X chromosome inactivation, and hematopoietic differentiation.

[Learn More](#)

Parameter	Value
Product Name	RBM15 Rabbit pAb
Remarks/Alias	OTT; OTT1; SPEN; RBM15
Species	Human
Gene ID (Human)	64783
Gene ID	64783
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 530-780 of human RBM15 (NP_001188474.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96T37
Protein Weight	107kDa
Shipping	Ice bag